DE PROTHESENKUNSTSTOFF

GB DENTURE BASE RESIN

**Xthetic[®]
cold**

**Xthetic[®]
cold
extra**

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Produktbeschreibung

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MADE IN GERMANY

Product description

Denture acrylic, autopolymer. Acrylic for fabricating removable or fixed dentures using the pouring technique.

Composition

2-component resin system consisting of powder and liquid. Main powder component: methyl methacrylate-copolymer. Main liquid component: methyl methacrylate, dimethacrylate.

Indications

- Fabrication of full dentures using the pouring technique
- Completion of partial dentures
- partial or full relining (indirect application)
- Additions and repairs

Contraindications

Patient contact with the unpolymerized acrylic or its individual components is contraindicated. Do not use in case of a known allergy to one or more components of the product.

Preparation

With partial dentures, the teeth set up in wax are fixed in position using a plaster or silicone index. Full dentures (wax setup on the model) are fixed in position in a suitable flask and the flask is filled with silicone or agar gel.

The teeth, metal framework and model are cleaned of wax residue by boiling and/or steaming and the model is soaked. With full dentures it is essential to soak the model until no more air bubbles rise from the plaster. This is particularly important to ensure that no air is pressed into the acrylic from the plaster during pressure polymerization which would result in porosities.

Separating

Coat the model and plaster indexes 1-2 times with a thin layer of alginate based separating agent. The **AcrylX® Xeparator** always shows reliable results.

Important: If separating twice, the first layer of separating agent should not be completely dry before applying the second layer. In certain circumstances this can result in white discoloration on the acrylic surface. It is essential to follow the instructions for use of the separating agent!

Note

Xthetic® cold and **Xthetic® cold extra** do not bond to highly cross-linked denture teeth. Denture teeth must be conditioned (roughened) and/or provided with mechanical retentions. It is recommended to coat the base surfaces of the teeth using **Xthetic® bond** to improve bonding strength.

Processing

Mixing ratio 10:7

The recommended mixing ratio is 10g powder and 7g liquid. Place the liquid in a mixing beaker, add the corresponding amount of powder as quickly as possible and mix thoroughly to a homogeneous consistency using a spatula.

Important: Thoroughly spatulate the mixture again after a saturation time of approx. 30 seconds to obtain homogeneous distribution of the saturated polymer beads! Ensure there are no bubbles in the mixture!

The prepared mixture is pourable for approx. 2-3 min. and can be filled into the indexes or openings of the casting canals.

With the beginning of the plastic phase, the resin becomes stable, no longer flows out of the indexes and can be modeled. The plastic phase lasts approx. 3 min for **Xthetic® cold** and approx. 20 min for **Xthetic® cold extra**.

Freehand dosing is possible, but the physical properties of the material are only guaranteed if the components are weighed out. When dosing freehand, add powder to the liquid until it produces a smooth consistency.

Polymerization - Xthetic® cold

Maximum 8 minutes after addition of the powder into the liquid, polymerization takes place under a pressure of 2-3 bar and a water temperature of 45-50°C [113-122°F] for approx. 7 min. The specification refers to a preceding processing time of approx. 8 minutes. The time in the pressure polymerization unit should last up to the 15th minute after mixing of the components.

Polymerization - Xthetic® cold extra

Maximum 20 minutes after addition of the powder into the liquid, polymerization takes place under a pressure of 2-3 bar and a water temperature of 45-50°C [113-122°F] for approx. 15-20 min. The specification refers to a preceding processing time of approx. 20 minutes. The time in the pressure polymerization unit should last up to the 35th minute after mixing the components.

After completion

Until incorporation the denture should be stored in water. It is always recommended to store the denture in water for 24 hours to achieve further reduction of the residual monomer content (2.4% according to ISO 20795-1), and thus prevent any taste irritation.

Feedback

Please direct any feedback regarding this product to **info@acrylx.com**. All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is located.

Storage and Safety Measures

Do not use the material after the expiration date. Store at temperatures not exceeding 25°C [77°F]. Protect from sunlight. Close containers immediately after use. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Dispose of contents/container in accordance with local disposal regulations.

Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/ face protection. IF ON SKIN: Wash with plenty of soap and water.

Colors

The material is available in colorless, pink, pink V (veined), pink TL (translucent), pink TLV (translucent-veined), pink opaque, V5 and pink C34.

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© = eingetragenes Warenzeichen der AcrylX[®] GmbH
® = registered trademark of AcrylX® GmbH